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© Meopant® Medical GmbH 2021

Vers. 01 – 17.Mai 2021

Instructions for use Meopant® Medical Implant system

English

Instructions for use

DISCLAIMER

This product is part of a comprehensive treatment concept and may only be used in combination with the associated original products in accordance with the instructions and recommendations of Meopant® Medical GmbH (Meopant). The non-recommended use of third party products in combination with Meopant products will void the warranty and void any other express or implied obligation of Meopant. The user of Meopant products must determine whether the product is suitable for a particular patient under the given conditions. Meopant accepts no liability, either express or implied, for any direct or indirect damages, punitive damages, or any other damages arising out of or in connection with any error in professional judgment or practice in connection with the use of Meopant products. The user is also obliged to regularly inform himself about the latest developments in relation to this Meopant product and its application. If in doubt, contact Meopant. As the use of the product is under the user's control, he assumes responsibility. Meopant accepts no liability whatsoever for damage resulting from the use of the product. We would like to draw your attention to the fact that some products mentioned in this manual may not be licensed or cleared for sale in all markets.

SCOPE OF THE INSTRUCTIONS FOR USE

These instructions for use apply to all components of the Meopant implant system, taking into account the Medical Devices Ordinance 2017/745.

The instructions for reprocessing the components that are approved for repeated use can be found in the instructions for use for preparation "Preparation instructions for the Meopant implant system" version 01 from 05/17/2021, which can be found at: <https://reprocessing.meopant.de/rp> and as a paper version with the delivered product(s).

safety instructions

Before using the Meopant® Medical implant system, the instructions for use for Meopant® Medical dental implants must be read carefully. The user is responsible for non-compliance with these instructions. Irrespective of all the information, the rules of dental surgery, occupational safety and accident prevention regulations must generally be observed for all the indications mentioned. Im-

proper use of the surgical and prosthetic parts can result in damage to the implant and the abutment parts, resulting in bone loss or even loss of the implant. Only properly trained dentists who are appropriately trained in implantology in relation to diagnosis, planning and surgical technology may use the Meoplast® Medical Implant System. The Meoplast® Medical Implant System is a coordinated implant system that must not be combined with structures and instruments from other implant manufacturers. In our terms of delivery, we guarantee the perfect condition and proper functioning of our products if all the above rules are observed and complied with. Before using the Meoplast® Medical implant system, the practitioner must ensure that all required product and system parts are in perfect condition and that they are secured against aspiration and swallowing during treatment. If there are any questions regarding the use or the indication limits before the operation, the operation must not be performed until these questions have been clarified. The practitioner concerned is solely responsible for any damage caused by the use of Meoplast® Medical products. The abutments/secondary parts may only be used according to their indication in accordance with the general rules for dental and surgical activities, taking into account the occupational health and safety and accident prevention regulations. Diameter reduced implants and angled abutments are not suitable for use in the posterior region. Advice on the use of Meoplast products takes place regularly and is strongly recommended.

DESCRIPTION

The implants from Meoplast® Medical are delivered sterile (gamma-sterilized) in double packaging. Before use, check the packaging for damage. The implants must be stored dry and at room temperature. If the expiry date has passed, these implants may no longer be used.

implant system


Figure 2: MeoClassic® and MeoMini® implant with closure screws

Meoplast® implants are made of grade 4 KV pure titanium. It consists of endosseous implants of different sizes (see Tables 1 and 2) and coordinated components, which are described below with regard to indication and application. The surface of the Meoplast® implants are blasted and acid-etched (SLA). The apically confiscated implant body has a cylindrical, parallel -walled shape with a macro-thread in the central part and a micro-thread in the crestal area. The three cutting edges ensure a self-centering insertion. The outer portions of the thread are designed so that a self-tapping effect occurs as long as the proposed drilling protocol is followed. The thread pitch and the furrow-like interruption of the thread has the function of collecting bone chips with vital bone cells, transporting them to the implant body and compacting them in order to biologically optimize the primary stability and healing into the bone. The tip of the implant has the shape of a lens and its convex area protects anatomically endangered structures such as the maxillary sinus mucosa in a sinus lift. The main purpose of the hexagon socket (hex connection) is to ensure a torsionally stable, permanently congruent connection between the implant and the abutment. The cone connection ensures an intimate metal connection, which reduces the micro gap (micro-gap) between abutment and implant and ensures better adhesion. The MeoClassic®

implant is designed in such a way that the implant edge should be flush with the crestal bone edge after insertion. The MeoTulip implant is designed in such a way that the polished tulip-shaped area should be visible above the crestal bone wheel after insertion. The diameter-reduced MeoMini® implant has an additional fanned, tulip-shaped polished shape instead of an internal hexagon, which merges into an external hex (Figure 2).


Figure 2: MeoMini • MeoClassic • MeoTulip - implant

The Meoplast® Medical Implant System consists of implants of different diameters and lengths intended for different indications. It is expressly pointed out that compliance with the prescribed drilling protocol is necessary. Separate components such as locking screws, gingiva formers, impression posts, MeoLock®, MeoMulti® and Meoplast® prosthetic components are available for the Meoplast® implant system. The MeoClassic and MeoTulip abutments / secondary parts are compatible with each other.

Abutments and secondary prosthetics

Figure 3: Abutments of the Meoplast implant system

Figure 4: Portfolio of screws of the Meoplast implant system

The Meoplast® Medical abutments can be used once and are supplied non-sterile. Before use, especially during the surgical phase of the dental implant procedure (the implant placement), they must be sterilized. Assemblies should be stored in a clean and dry area.

instruments


Figure 5: Portfolio of instruments of the Meoplast implant system


Figure 6: Tools for taking an impression of the Meoplast implant system

PURPOSE OF USE

The Meoplast implants are dental implants that are intended for use in the upper or lower jaw bone and are used to anchor or attach dentures to restore chewing function.

The drills and accessories of the implant system are instruments for the application of endosseous Meoplast dental implants, which are used in the event of the loss of one or more teeth to restore the human dentition by means of a fixed denture.

INDICATIONS

The Meoplast® Medical implant system is intended solely for use in dental implantology and maxillofacial surgery. For single crowns and bridges, both screwed and cemented, to support partial dentures, full dentures and epitheses. In all applications, taking into account the bone quality, always aim for axial loading and a sufficient implant length, which should not be less than 1:1 in relation to the crown length. The Meoplast implant system is intended for surgical procedures in the sense of immediate, delayed and late implantations, both single-phase and two-phase, for interstitial gaps, shortened rows of teeth and edentulous jaws. The diameter-reduced MeoMini®, Ø2.9mm implant, is not suitable for the single-tooth restoration of central incisors in the upper jaw and canini, premolars and molars in the lower and upper jaw. In this region (upper jaw) the implant sizes Ø4.2mm, Ø4.8mm and Ø6.0mm must be used. In the lower jaw, the Ø3.8mm can also be used in the side teeth area. However, the implant length should not be less than 10mm or, in the case of shorter implants, should be splinted. The MeoMini® implant is suitable for increasing pillars and for supporting prostheses and bridges, as well as for use in orthodontics as a counter bearing.

CONTRAINDICATIONS

Wound healing disorders, endocrine diseases, osteomyelitis, systemic diseases of the bones, the immune system, the haematopoietic system, rheumatoid diseases, liver diseases, diabetes that is difficult or impossible to control, mental diseases. Temporary contraindications: nicotine and alcohol abuse, poor oral hygiene, acute inflammation in the surgical area, insufficient bone volume and endangerment of important anatomical structures, recurrent oral mucosal diseases, acute parafunctions, lack of patient compliance, insufficient interocclusal distance, insufficient soft tissue coverage, steroid therapy, Radiation therapy in the head area, chemotherapy, drugs that interfere with the calcium balance, anticoagulant therapy.

Products that are not part of the Meoplast® Medical Dental implant system must not be used.

SIDE EFFECTS AND COMPLICATIONS

The usual side effects of surgical interventions such as swelling, hematoma, edema and postoperative pain, restricted mouth opening and chewing as well as possible sensory disturbances are possible.

In rare cases, the following complications can occur intraoperatively:

- infections of the surgical area,
- suture dehiscence,
- insufficient primary stability,
- unwanted trauma,
- postoperative bleeding from vascular injuries,
- injuries to important anatomical structures,
- broken screws,
- canting of locking screws or gingiva formers,
- Break or fracture of the implant neck and
- Aspiration or ingestion of small parts used.

Allergic reactions to titanium are extremely rare, but cannot be ruled out. In special cases, reactions of different alloy components of the implant and abutments can lead to allergic reactions. This can cause galvanic reactions within the oral cavity, which can affect the patient's well-being. The practitioner must ensure that complications of this type are avoided.

General

There is no 100% guarantee of success for implants. Failure to comply with the specified limitations of use and work steps can lead to functional errors.

The insertion of implants can lead to bone loss and biological or mechanical failure, eg fatigue fracture of the implant.

The close cooperation between surgeon, prosthodontist and dental technician is essential for successful implant treatment.

Improper surgical and prosthetic procedures can result in damage to the dental implant or loss of bone. The Meoplast implant system may only be used by dentists, doctors and surgeons trained with the system, since the use of the Meoplast implant system requires special knowledge and skills in implantology.

The components of the Meoplant implant system are only given to dentists/MCG and dental laboratories or on their behalf. This is to ensure that the specific knowledge that enables safe application is available.

Inhomogeneities can be visually observed on the surface of the implant. Complete wetting of the surface is ensured microscopically. Any inhomogeneities that occur therefore have no impact on product quality, clinical safety and clinical effectiveness.

Only the tools of Meoplant Medical GmbH may be used for the insertion. A torque wrench is available for this purpose. This can be equipped with all the necessary tools via the ratchet adapter. If you are working with a contra-angle that has an ISO shank, the tools can also be used with these contra-angles. The use of non-system components and instruments can impair the function and safety of the Meoplant implant system. If non-system components are used, no guarantee and no replacement will be provided.

Drills, instruments and system components are designed for specific implants and implant diameters. The different diameters are identified by color markings (see Tables 1 and 2 and Figure 7). Use for other implants or other diameters may result in mechanical failure of system components, tissue damage, or unsatisfactory esthetic outcomes.

Table 2: Coding of the MeoClassic implant system

implant	Ø [mm]	length [mm]	colour
MeoMini®	2.9	8 / 10 / 12.5	White
MeoClassic®	3.5	8 / 10 / 12.5 / 15	yellow
MeoClassic®	3.8	8 / 10 / 12.5 / 15	red
MeoClassic®	4.2	8 / 10 / 12.5 / 15	green
MeoClassic®	4.8	8 / 10 / 12.5 / 15	blue
MeoClassic®	6.0	8/10	black

Table 2: Coding of the MeoTulip (Tissue) implant system

implant	Ø [mm]	length [mm]	colour
MeoTulip®	3.5	8 / 10 / 12.5	yellow
MeoTulip®	3.8	8 / 10 / 12.5	red
MeoTulip®	4.2	8 / 10 / 12.5	green
MeoTulip®	4.8	6 / 8 / 10 / 12.5	blue
MeoTulip®	6.0	6 / 8 / 10	black


Figure 7: Recommended implant diameters in relation to tooth position

Before the intervention

Before the surgical implantation, a careful patient history and, if necessary, a third-party history must be taken by the family doctor to clarify whether

- e) a more difficult implantation due to the anatomical situation,
- f) a serious surgical problem or general risk,
- g) impaired wound healing and/or osseointegration, or
- h) deterioration in proper hygiene and/or care of the implant, abutment and prosthesis could occur.

During the procedure

All instruments and equipment used during the procedure must be in good condition and care must be taken to ensure that the instruments do not damage the implants or other components.

Due to the small size, choking and aspiration of a component of the Meoplant implant system can occur. Aspiration can lead to shortness of breath and, in the worst case, to asphyxiation. For this reason, the products must be secured against swallowing and aspiration with a thread or dental floss when used intraorally.

It is essential to follow the drilling sequence while preparing the implant bed. The recommended torques are described in Table 3.

Table 3: Recommended torques of the superstructures / secondary components

system part (superstructures / secondary parts)	Torque [Ncm]
Meoplant® straight & angled / angled	30
Meoplant® one-piece solid / ball / bar system	
MeoLock® one-piece & two-piece	
Perform / Cerec ® Adhesive Base	
MeoMulti® straight & angled / angled	27
adhesive base	
MeoLock® Cap	
MeoMulti® structure (titanium)	
MeoMulti® safety cap	25
MeoMini® one-piece & two-piece / ball head	
PEEK & PEKK (High Performance Thermoplastics)	15
Impression posts / POM abutments	5-10
Screw plug / cover screw	2-4
Gingiva formers / healing caps	

After inserting the implant, the surgeon's assessment of bone quality and primary stability determines when the implant can be loaded. Inadequate quantity and/or quality of available bone, infections and systemic diseases are among the possible reasons for a lack of osseointegration, both immediately after the procedure and after an initial osseointegration.

If the implant or a screw-in head fractures during insertion due to excessive stress on the material, both can be removed again using the rescue tool. To do this, attach the rescue tool of the appropriate size and remove the implant and the screw-in head with a left-turning movement.

After the procedure

To ensure optimal long-term treatment outcome, regular comprehensive follow-up appointments should be scheduled with the patient and the patient should be informed about optimal oral hygiene.

SURGICAL PROCEDURE

Care must be taken during all steps not to compromise the sterility of the product and the sterile or sterilized product components.

Every implantation is a surgical procedure that requires education, planning, sterile precautions and gentle treatment of the soft and hard tissues. This includes drilling with sufficient cooling and with

no more than 400 rpm,

with intermittent movement and low pressure. Multiple use of drills should be avoided, as there is a risk of bone overheating if drills are too blunt and healing (osseointegration) in the bone cannot occur due to heat necrosis. After preparing the bone cavity, the MeoClassic (1) / MeoTulip (1.1) implant is inserted with the insertion tool (56 | 57) and the MeoMini® implant (47) with the insertion tool (63) at 10-15 rpm. The insertion torque when inserting the implant should not exceed 50 Ncm and under no circumstances should it be higher in order to avoid pressure necrosis and thermal damage. Likewise, make sure that the screw-in torque does not fall below 15 Ncm. In general, the insertion tools should always have an exact fit to the inner geometry of the implant when being used in order to avoid damaging the inner geometry of the implant. The implant surface is specially produced as an SLA-compatible surface and is highly active. All touching of the implant surface prior to insertion should be avoided. The MeoClassic® implant (1) is constructed in such a way that the implant edge should be flush with the crestal bone edge after insertion. The MeoTulip implant (1.1) is designed in such a way that the polished tulip-shaped area should be visible above the crestal bone wheel after insertion. In addition, the exact prosthetic position of the subsequent crown should match the position of the implant in the horizontal and sagittal direction through precise operation planning, taking into account the bone availability. An implantation can be carried out with the formation of a mucoperiosteal flap or also transgingivally, this requires precise knowledge of the bone volume. In the case of a transgingival procedure, the soft tissue of the implant must be avoided at all costs. The healing can be closed after screwing in a locking cap/screw (2) and suturing the incision slot or in an open way after screwing in a healing cap (3 | 4) or a ready-made, possibly previously individualized abutment with a long-term temporary made for it in a perforated with a gingiva punch sulcus or a flap sewn all around. There must be neither an approximal nor an occlusal contact with neighboring teeth.

Before inserting the implant, the depth gauge can provide orientation about the drill channel depth. The actual depth can only be determined by selecting the drills and the implants. The depth gauge is only for orientation during the surgical procedure.

Exposure / impression / prosthetic care

In the case of single-phase healing/implantation, the implant (1 | 1.1) is fitted with a healing cap/gingiva former (3 | 4) or an abutment (10 - 16 | 17 - 22 | 27 | 33 | 34) after insertion. Immediate restoration with a temporary restoration on a prepared abutment after an impression can also be made. It is important to ensure that there is no approximal or occlusal contact with neighboring teeth. In the two-phase procedure, the implant closure cap/screw (2) is exposed and removed after the specified healing time under anesthesia. After optionally closed (6 / 7) or open impression (5), gingiva formers (3 | 4) or the respective abutment are reinserted. Impression posts (5 | 6) must be screwed in with a maximum of 5-10 Ncm, whereby the exact fit is a prerequisite for the fitting accuracy of the prosthetic restoration. For open impression taking (5), a passage opening for the impression post (5) is ground into the individual impression tray so that the impression post and implant can be detached from one another after the impression has been taken. For further processing in the laboratory, the impression post (5) in the impression is screwed to a laboratory analog (8 or 9) that has the same inner geometry as the implant. After the model has been made, the screw connection of the laboratory analogue and the impression post is separated. With the closed impression (6 + 7), the corresponding impression post (7) is connected to the implant body (1 | 1.1). A transfer cap (6) is then placed on the impression post in the mouth. After the impression material has set and the tray has been removed, the transfer cap (6) is located in the impression material of the tray. The impression post (7) is now detached from the implant and connected to a laboratory analogue (8 | 9). The model is made in the same way as for open impressions. After the individualized abutment has been made in the laboratory and, if necessary, the definitive denture, the provisional integration takes place. The peri-implant tissue maturation takes about six months and is used for the long-term prognosis of the implant. A laboratory screw (50) (blue anodized) specially developed for this work process must be used for all work in the laboratory. When inserting the final prosthetic work, use the definitive fixation screw (51) provided with the relative abutment, screwing either just the abutment or the whole crown. Cementing is an alternative to screwing. When selecting the abutment, attention must be paid to the height of the gingiva, the diameter and the angulation of the implant. A fracture of the connecting screw can be avoided by observing the maximum torque of 30 Ncm. However, should a screw break, the fractured screw is removed with a specially designed rescue kit. The adhesive labels included in the

packaging with the specific article information are available for documentation.

closed healing

implant insertion

Before the suture, the internal geometry is closed with the help of a locking screw (2 | 45), which is screwed in hand-tight with a maximum of 2 - 4 Ncm with the hexagonal instrument (58 | 59 / 60 | 61), since screwing in too tightly will result in later loosening of the screw plug during exposure is made more difficult or impossible.

exposure

With the hexagonal instrument (1.2 58 | 59) for MeoMini® (47) or with the hexagonal instrument (1. 4 60 | 61) for the implants (1|1.1), the cover screw (2 | 45) is removed after exposure and provide the implant with a healing cap (3 | 4). However, an impression can be taken before or after the soft tissue conditioning phase.

b) open impression

Impression body (5) with fastening screw (55) is screwed onto the implant (1 | 1.1) with a maximum of 5-10 Ncm in a short or long version, depending on the space available and localization. It is important to ensure that there are no tissue or foreign particles in the inner geometry of the implant in order to ensure the exact force and form fit. A passage opening is provided for the fastening screw in a previously prepared individual impression tray in order to loosen the screw from the implant (1 | 1.1) after the impression material has set and to be able to pull it out of the impression body (5).

b) closed impression

A three-part impression post (7 | 37) is available for closed impressions. The impression body is screwed to the implant (1 | 1.1) with the aid of the fastening screw (51) with a maximum of 5-10 Ncm and fitted with the corresponding cap (6), which audibly clicks into place on the impression body. Conventional impressions can now be taken with a closed tray, although alginate is generally not suitable for implant impressions. The impression cap 6 remains in the set impression compound and, after it has been released, enables the impression body connected to a model analog (8 | 9 | 38) to be repositioned in the impression.

Laboratory - model making

The technician again checks the exact fit of the respective impression body in the impression material before selecting the appropriate laboratory analogue (8 | 9) and fastening it to the impression post (5 | 7) with the laboratory screw (50). After the application of a gingiva mask made of hardening rubber and the actual model production with plaster of paris, the intended abutments can be further processed.

Laboratory – abutment selection and processing

Depending on your wishes or the axis of the placed implant, abutments (10 - 16 | 18 - 22 | 27 | 33 | 34 | 46 | 48 | 49) can be selected and individualized from the selection offered. Standard abutments (10 | 12) are sufficient if the implant is aligned correctly. For minor corrections, the one-piece structure (13) can be adjusted to the tooth axes by grinding in the mouth using a rubber dam. The universal abutments (15 | 16) can be individualized in the laboratory.

open healing

suitable healing caps made of titanium (3) with 2-4 Ncm or customizable healing caps made of Peek (4) with 15 Ncm can be screwed onto the connecting screw (5) with open healing, depending on the thickness of the gingiva, in order to close the peri-implant soft tissue to form and mature. Impression as described under a) and b).

Open healing with abutments and long-term temporary

If the soft tissue maturation is to proceed as undisturbed as possible, frequent changing of structural parts is contraindicated. Here it makes sense to insert the abutments in the same way as the later dentures. Depending on the fit of the implant, straight abutments (10 | 11 | 13 - 16) can be used, which can be individually adjusted to the tissue conditions by grinding .

Open healing with abutments and prefabricated dentures for immediate loading

The only difference from the previous paragraph is the supply. As part of CAD/CAM planning, the dental technician can then digitally scan the dentures using scan abutments (17 | 43), which are screwed onto the implant that has just been placed in the mouth. With the resulting data set, the crown, bridge or bar restoration

can be designed and used accordingly in the dental technology laboratory. At least four implants are required for immediate loading.

MeoLock®

Denture prosthesis is secured against vertical pull-off forces by the MeoLock® abutment system . One-piece and two-piece MeoLock® abutments (21 + 22) are available in heights of 1-6 mm. After exposure, the one-piece MeoLock® abutment (21) is screwed onto the implants (1 | 1.1) with the MeoLock® insertion tool with 30 Ncm. The two-part MeoLock® abutment (22) is screwed onto the implant (1 | 1.1) with the connecting screw (51) with 30 Ncm and provided with an impression cap (23) for the impression. After the impression has been taken, the dental technician places the MeoLock® model analog (24) in the impression cap (23) remaining in the impression and fabricates the master model. The cover denture prosthesis is now constructed over the MeoLock® attachment (21 | 22). The MeoLock® matrix (25) can be polymerized both in the mouth and on the model . The colors of the retention elements (26) (patrix) symbolize the different pull-off forces. If the retention elements are worn , they can be replaced.

ball head attachments

Ball head abutments (27) with the Insertion instrument (62) with 30 Ncm and with the MeoMini® implant (47) with 25 Ncm. A ball head model analogue (28) is available for model production. For the impression, the ball head abutments (46) for MeoMini® (47) and the ball head abutments (27) for MeoClassic (1) and MeoTulip (1.1) are screwed on and conventionally molded. After removing the impression tray, the ball head analogue (28) is inserted into the resulting cavity and the master model is made in the laboratory. With both types of implants , the dental technician has the option of individually fitting the two existing ball head abutments (27 | 46) with two different retention abutments. Firstly, the rubber O-ring (30) and the associated O-ring mount (29) as a matrix that is polymerized into the prosthesis and secondly, the Preci - Clix male (29) with the associated Preci - Clix matrix (31) which is also polymerized into the prosthesis .

MeoMini® bodies

The one-piece MeoMini® compact abutment (48) is screwed onto the MeoMini® implant (47) in the mouth with the insertion instrument (60 | 61) with 25 Ncm and then molded. The two-part MeoMini abutment (49) is fixed to the MeoMini® implant (47) with the connecting screw (53) with 25Ncm and then molded.

bar system

Three versions are available for the laboratory production of bar constructions. The gold bar cylinder (18) is used by the dental technician for casting. The plastic bar cylinder (19) is embedded together with the modeled bar construction and cast in metal. This plastic cylinder can be burned out. Another bar cylinder is made of titanium (20) and can be fixed with a laser in the laboratory. The selected bar abutments are fixed on the model with the laboratory screw (50) and can now be processed. After completion of the bar construction, it is fastened in the mouth with the final connecting screw (51) with 30Ncm.

CLEANING AND STERILIZING INSTRUCTIONS

All MeoPlant implants are delivered sterile and are intended for single use only before the stated expiry date.

Warning: The components that are delivered sterile must no longer be used if the packaging is damaged or has been previously opened.

INFORMATION ON MRI SAFETY

In the case of magnetic resonance examinations, it should be noted that the MeoPlant implant system has not been tested for safety and compatibility in MRI examinations.

HANDLING AND STORAGE

The implants are gamma sterilized in the packaged condition and are intended for single use. There must be no resterilization. The implants may no longer be used if the packaging is damaged or opened and if the expiry date has expired. The manufacturer declines all responsibility in the event of non-compliance. All instruments and abutments/secondary parts are packaged non-sterile and must be sterilized before each use ("Preparation instructions MeoPlant implant system" version 01 from May 17th, 2021 can be found at: <https://www.meoPlant.info/download/sa/mobile/index.html> and as a paper version with the delivered product). For all items intended for single use, it was not tested whether repeated cleaning procedures and sterilization processes have an effect on the function and fitting accuracy of the material. Therefore, multiple use of this item is strongly discouraged. The product must be stored in a dry place in its original packaging, at room temperature and away from direct sunlight. Improper storage can affect product properties and lead to supply failure.

DISPOSAL

The product must be disposed of in accordance with local regulations and environmental regulations, taking into account the degree of contamination.

Manufacturer:

Meoplant Medical GmbH
Malchiner Strasse 99 – 12359 Berlin / Germany
Tel. +49 30 809334166 / info@meoplant.de /
<https://www.meoplant.de>

EXPLANATION OF SYMBOLS

Do not resterilize

Do not reuse

Attention: Observe the accompanying documents

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- All rights reserved -
© Meoplant® Medical GmbH 2021